

PRIMA C/CS



A POWERFUL SOLUTION IN A COMPACT DESIGN




The PRIMA C is a compact and versatile solution for surgical and diagnostic work in Gynecology. Bundled with a high power LED illumination module, superb optics and an ergonomic footprint, this colposcope ensures durability and high performance.

The PRIMA C comes mounted on a choice of superbly crafted stands, making it adaptable to any space or application. A retractable pole mount with an elegant star base, a foldable suspension arm, and various chair mounts make PRIMA C systems highly customizable.

A range of LABOMED accessories *Proline* cameras & smart phone adapters enable image streaming and archiving for any Gynecology practice.

A high performance LED illumination system built into the magnichanger provides superb lighting of the area under observation. LED light levels exceed those provided by halogen lamps, making the PRIMA C design compact and well balanced. This unique arrangement is fitted with an ultra bright LED rated for up to 50,000 hours, virtually eliminating the need for bulb changes.



LABOMED's proprietary *MaxLite™* coated apochromatic optics deliver a level of detail, color correction, and anti-reflective properties typically found in operating room microscopes, making the PRIMA C a best-in-class value colposcope.

PRIMA C COLPOSCOPE

Observation tubes

IPD

Eyepieces

Dioptre Adjustment

Magnichanger

Objective

Magnification

FOV

Light Source

LUX

Maximum Wattage

Color Temperature

Built-in Filters

Vertical Movement of Arm

Optional accessories

PRIMA C Colposcope with Straight Carrier

Standard: Binocular straight 90° viewing tube
Optional: 45° Binocular inclined

49-78 mm

Standard: WF 10X/18 mm with eye guards

Option: WF 12.5X/18 mm with eye guards

10x: ± 7 mm, 12.5x: ± 8 mm

Apochromatic *MaxLite* coated

5 step: 0.4X, 0.6X, 1.0X, 1.6X, 2.5X

Standard: f-300 mm, manual focus

Option: f-400 mm, manual focus

2.3 x-14x

79.5 - 12.70mm

27W LED, 100-240V

75K

60 Watts

5700k

Green

100 mm

Beam splitter (straight), Double Iris,

ProLine range of camera adapters &

Smart phone adapters,

Assistant binocular attachment

PRIMA CS Colposcope with Swing Arm

Standard: Binocular Straight 90° viewing tube
Optional: 45° Binocular inclined

49-78 mm

Standard: WF 10X/18 mm with eye guards

Option: WF 12.5X/18 mm with eye guards

10x: ± 7 mm, 12.5x: ± 8 mm

Apochromatic *MaxLite* coated

5 step: 0.4X, 0.6X, 1.0X, 1.6X, 2.5X

Standard: f-300 mm, manual focus

Option: f-400 mm, manual focus

2.3 x-14x

79.5 - 12.70mm

27W LED, 100-240V

75K

60 Watts

5700k

Green


± 250 mm

Beam splitter (straight), Double Iris,

ProLine range of camera adapters &

Smart phone adapters,

Assistant binocular attachment



ITEM

PRIMA C

CONFIGURATION

Colposcope with
Straight Carrier

CATALOG NO.

6128000

PRIMA CS

Colposcope with
Swing Arm

6128600



LABOMED®

Labotech Microscopes India Pvt Ltd
#19 HSIDC Industrial Estate
Ambala Cantt
Haryana, Pincode: 133006
INDIA

Distributor