

ERGOFLEX




PERFECT COMFORT - PAIN FREE DAY

Ergoflex is the result of years of working closely with doctors and understanding causes of fatigue and lumbar stress due to prolonged hours of working in uncomfortable postures. The way you sit is no minor thing – pain is what an uncomfortable posture will bring.



SPECIALY ENGINEERED SUPPORTIVE BACKREST

The special spring loaded backrest provides constant lumbar support greatly reducing fatigue to the doctor's back in every sitting posture. Data supports that 8 out of 10 dentists suffer from pain caused by long hours of uncomfortable seating. Ergoflex presents a solution that is good form, good function translating to good health.




OPTIMAL ADJUSTABILITY

Independently adjustable armrest provides comfort customisable to every working style and body type. Foot controlled height adjustment provides hands free ergonomic set up. Your chair is your work zone – adjust it to fit your needs alone.




FREEDOM OF MOVEMENT FREEDOM FROM DISCOMFORT


Loaded with unmatched features, Ergoflex provides a quick and easy seating solution for doctors and their assistant. Adjustable height, seat angle and elbow rest settings features prevent Musculoskeletal disorder. A seating system to ensure that you are not bent out of shape.



Ergo flex is also particularly hygienic thanks to the possibility it offers with hands free seat height adjustment by additional foot control mechanism. Upholstery covering of Ergoflex is easy to clean and anti-microbial.



Five leg spider base with double castor wheels for high stability, safe, precise and easiest movement on all hard floors.




Ergoflex impresses doctors and assistants alike. With reliable quality, modern design, excellent price-performance ratio, Ergoflex makes for an excellent choice for all operatories.



Grey Maroon Navy Blue Black

Available in a choice a colours to complement various decors and applications.

GOOD FORM. GOOD FUNCTION. GOOD HEALTH.



ErgoFlex

Seat height	Minimum 480mm
Extension	Upto 660mm
Back rest tilt	0° - 25°
Back rest height adjustment	250mm
Arm rest height adjustment	150mm
Arm rest width	Minimum 530mm
Extension	Upto 735mm
Arm rest radial adjustment	± 50°
Base diameter	560mm
Fabric colour	Black, Grey, Maroon, Navy Blue
Weight	24 kg
Seat system design	Contoured water fall design with 110 degree angle
Pneumatic cylinder	Anticorrosive steel, 150 kg capacity
Upholstery	High resilience PU foam, anti fading & anti UV fabric
Base and arm rest	Made from high polished anti rust Stainless Steel

SPECIFICATIONS



A16372 ISO 13485 : 2016



LABOTECH MICROSCOPES INDIA Pvt Ltd

Labotech Microscopes India Pvt Ltd
#19 HSIDC Industrial Estate
Ambala Cantt
Haryana, Pincode: 133006
INDIA



6137160-796/E/ 1.0 010822